

# ● Analysis Report

## Bruker IVDr Quantification in URine B.I.Quant-UR b<sup>TM</sup>

Sample ID: ALZ\_Urine\_Rack01\_RCM\_221214\_expno440.100000.10r

Measuring Date: 24-Dec-2014 00:18:19

Reporting Date: 12-Dec-2020 13:27:09, 7 page(s), Version 1.1.0

Quantification Method Version: Quant-UR B.1.1.0

### Disclaimer

RESEARCH USE ONLY: This is no clinical diagnostic analysis report. Must not be used for clinical (medical or IVD) diagnosis or for patient management! Additional concentration range information (95% range) provided numerically or graphically in this report must not be used for clinical diagnostic interpretation.

Application of B.I.Quant-UR B 1.1.0 requires use of Bruker's B.I.Methods SOP for urine.

### Summary

The following metabolites were found with concentrations outside the 95% range of Bruker Quant-UR B.1.1.0 urine metabolite concentration database:

Benzene and substituted derivatives: Hippuric acid (940 mmol/mol Crea),


Carboxylic acids: Acetic acid (220 mmol/mol Crea), Fumaric acid (4 mmol/mol Crea), Succinic acid (44 mmol/mol Crea), Tartaric acid (200 mmol/mol Crea), Trigonelline (85 mmol/mol Crea).

Further detailed information is provided on the following pages.

## Contents



<b>1 Creatinine</b>	<b>3</b>
<b>2 Amines and derivatives</b>	<b>3</b>
<b>3 Amino acids and derivatives</b>	<b>3</b>
<b>4 Benzene and substituted derivatives</b>	<b>4</b>
<b>5 Carboxylic acids</b>	<b>4</b>
<b>6 Fatty acids and derivatives</b>	<b>4</b>
<b>7 Keto acids and derivatives</b>	<b>5</b>
<b>8 Purine, Pyridine and Pyrimidine derivatives</b>	<b>5</b>
<b>9 Sugars and derivatives</b>	<b>5</b>
<b>10 Explanations</b>	<b>6</b>

## 1 Creatinine

Compound	Conc. mmol/L	LOD mmol/L	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range <sup>(*)</sup> mmol/L
Creatinine	3.0	0.3	2.973	100 ●	0.069	1 - 19 















(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 2 Amines and derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range <sup>(*)</sup> $\frac{\text{mmol}}{\text{mol Crea}}$
Dimethylamine	0.09	31	31	0.092	100 ●	0.001	$\leq 54$ 
Trimethylamine	< 0.01	< 2	2	0.003	52 ○	0.002	$\leq 3$ 




(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 3 Amino acids and derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range <sup>(*)</sup> $\frac{\text{mmol}}{\text{mol Crea}}$
1-Methylhistidine	< 0.04	< 15	15	0.000	0 ○	0.147	$\leq 15$ 
2-Furoylglycine	< 0.11	< 39	39	0.026	80 ○	0.008	$\leq 40$ 
4-Aminobutyric acid	< 0.06	< 20	20	0.022	35 ○	0.029	$\leq 20$ 
Alanine	0.06	20	10	0.060	100 ●	0.004	11 - 72 
Arginine	< 2.2	< 750	750	0.167	0 ○	0.648	$\leq 750$ 
Betaine	0.04	15	7	0.043	100 ●	0.006	9 - 78 
Creatine	< 0.15	< 50	50	0.068	100 ●	0.069	$\leq 280$ 
Glycine	0.25	83	34	0.248	100 ●	0.005	38 - 440 
Guanidinoacetic acid	< 0.31	< 100	100	0.143	86 ●	0.047	$\leq 140$ 
Methionine	< 0.05	< 18	18	0.000	0 ○	0.080	$\leq 18$ 
N,N-Dimethylglycine	0.02	6	5	0.018	77 ○	0.009	$\leq 15$ 
Sarcosine	< 0.01	< 2	2	0.003	0 ○	0.005	$\leq 7$ 
Taurine	< 0.42	< 140	140	0.259	59 ○	0.131	$\leq 170$ 
Valine	< 0.01	< 2	2	0.000	0 ○	0.039	$\leq 7$ 











(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 4 Benzene and substituted derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
Benzoic acid	< 0.03	< 10	10	0.000	0○	0.031	≤ 10 
D-Mandelic acid	< 0.01	< 2	2	0.000	0○	0.057	2 - 17 
Hippuric acid	2.8	940	170	2.791	99●	0.270	≤ 660 


(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 5 Carboxylic acids

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
Acetic acid	0.65	220	5	0.653	100●	0.012	≤ 51 
Citric acid	0.70	240	40	0.703	100●	0.067	≤ 700 
Formic acid	< 0.03	< 10	10	0.012	93●	0.003	≤ 43 
Fumaric acid	0.01	4	2	0.012	100●	0.000	≤ 3 
Imidazole	< 0.14	< 48	48	0.000	0○	0.070	≤ 48 
Lactic acid	< 0.14	< 49	49	0.064	46○	0.066	≤ 110 
Proline betaine	< 0.07	< 25	25	0.048	0○	0.087	≤ 280 
Succinic acid	0.13	44	5	0.131	100●	0.008	≤ 39 
Tartaric acid	0.60	200	5	0.602	100●	0.014	≤ 110 
Trigonelline	0.25	85	35	0.252	100●	0.006	≤ 67 







(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 6 Fatty acids and derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
2-Methylsuccinic acid	< 0.14	< 48	48	0.000	0○	0.153	≤ 48 








(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 7 Keto acids and derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
2-Oxoglutaric acid	< 0.27	< 92	92	0.029	36 ○	0.074	≤ 92 
3-Hydroxybutyric acid	< 0.31	< 100	100	0.000	0 ○	0.195	≤ 100 
Acetoacetic acid	< 0.04	< 14	14	0.029	81 ○	0.022	≤ 30 
Acetone	0.01	3	2	0.009	99 ●	0.001	≤ 7 
Oxaloacetic acid	0.06	20	17	0.060	82 ○	0.043	≤ 66 
Pyruvic acid	< 0.03	< 9	9	0.014	92 ●	0.006	≤ 13 







(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 8 Purine, Pyridine and Pyrimidine derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
1-Methyladenosine	< 0.01	< 5	5	0.000	0 ○	0.065	≤ 5 
1-Methylnicotinamide	< 0.09	< 32	32	0.049	98 ●	0.007	≤ 32 
Adenosine	< 1.2	< 390	390	0.000	0 ○	0.701	≤ 390 
Allantoin	< 0.05	< 17	17	0.021	93 ●	0.005	≤ 47 
Allopurinol	< 0.03	< 10	10	0.020	92 ●	0.024	≤ 11 
Caffeine	< 0.13	< 45	45	0.083	80 ○	0.075	≤ 61 
Inosine	< 0.06	< 19	19	0.012	96 ●	0.027	≤ 19 

(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 9 Sugars and derivatives

Compound	Conc. mmol/L	Conc. $\frac{\text{mmol}}{\text{mol Crea}}$	LOD $\frac{\text{mmol}}{\text{mol Crea}}$	r mmol/L	$\rho$ %	$\Delta$ mmol/L	95% Range(*) $\frac{\text{mmol}}{\text{mol Crea}}$
D-Galactose	< 0.13	< 43	43	0.000	0 ○	0.014	≤ 44 
D-Glucose	< 0.10	< 34	34	0.029	58 ○	0.021	≤ 140 
D-Lactose	< 0.28	< 96	96	0.154	89 ●	0.062	≤ 96 
D-Mannitol	< 0.54	< 180	180	0.000	0 ○	1.270	≤ 180 
D-Mannose	< 0.02	< 6	6	0.007	0 ○	0.005	≤ 8 
Myo-Inositol	< 13	< 4400	4400	0.000	0 ○	1.603	≤ 4400 

(\*) Gray horizontal boxes represent 95% concentration range, black vertical lines represent sample value.

## 10 Explanations

This section contains the definition of the parameters used above. In the section 10.1 a short manual, how to interpret the results, is presented. The section 10.3 contains the exact definitions of the parameters  $r$ ,  $\rho$  and  $\Delta$ .

### 10.1 How to read the result



Figure 1: Examples of fitting.

In the figure 1(a), the black line, the blue line and the yellow line represent the original spectrum, the calculated signal fit and its baseline, respectively.

The blue area relates to the metabolite concentration to be determined and the red area represents a residue.

In case of the signal overlap a different approach is used: two or more overlapping signals are being fitted simultaneously. The most iconic example of such signals are the ones generated by  $\text{CH}_3$  groups of Creatinine and Creatine. In such a case, the blue line and the grey area relate the sum of all fitted signals. The blue area corresponds to the concentration of the metabolite of interest (cf. figure 1(b)).

### 10.2 Result parameters

- Conc.** is the final result concentration of the metabolite,
- LOD** is the *limit of detection* of the given metabolite,
- r** is the *raw concentration* i.e. the concentration equivalent of the resulting signal fit prior to comparing to **LOD** (relates to the blue area, cf.  $\alpha$ ),
- $\rho$  is the correlation of lineshape metabolite signal with calculated fit characterizing the match between metabolite signal and fit (cf.  $\beta$ ). Depending on the value of  $\rho$ , the following *flag* is displayed:

- ●, if the correlation is 95%,
  - ●, if the correlation is in between 85% and 95%,
  - ○, if the correlation is less than 85%,
- e)  $\Delta$  is the concentration equivalent of the difference between metabolite signal and calculated fit (residue corresponding to the **the red area**, cf.  $\gamma$ )).

### 10.3 Detailed definitions

Let  $s$ ,  $f$  and  $b$  denote the functions describing the *raw spectra*, *fitted curve* and *(fitted) baseline* respectively. These functions are chosen such that  $s \approx f + b$ . Moreover, let  $I$  be a relevant PPM interval and  $P_N$  be the proton number for given metabolite/signal.

$\alpha$ )  $r$  (*raw concentration*) is defined as

$$r = \frac{1}{P_N} \int_{\mathbb{R}} f(\xi) d\xi.$$

$\beta$ )  $\rho$  is the *correlation* of the functions  $s$  and  $f + b$ , i.e.

$$\rho = \max(0, \text{corr}(\bar{s}, \overline{f+b})),$$

where  $\bar{s}$ ,  $\overline{f+b}$  are numerical representations of the functions  $s$  and  $f + b$  on sufficiently fine mesh of the interval  $I$ .

$\gamma$ )  $\Delta$  is the the area between the raw signal  $s$  and the fitted data  $f + b$  on the interval  $I$  expressed in the terms of the concentration, i.e.

$$\Delta = \frac{1}{P_N} \int_I |s(\xi) - f(\xi) - b(\xi)| d\xi.$$